
## Single cell RNA Sequencing 분석 기술

### < Single Cell RNA-Seq vs Bulk-Seq >

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### Single cell RNA-Seq

### Bulk RNA-Seq

<b>Goal</b>	<ul style="list-style-type: none"> <li>• Measure the gene expression of individual cells in a sample</li> <li>• To identify differences between <b>cell types/states</b></li> </ul>	<ul style="list-style-type: none"> <li>• Measure the average gene expression across the population of cells in a sample</li> <li>• To identify differences between <b>sample conditions</b></li> </ul>
<b>Feature</b>	<b>Separate Populations</b> <ul style="list-style-type: none"> <li>• Define heterogeneity</li> <li>• Identify rare cell populations</li> <li>• Cell population dynamics</li> </ul>	<b>Average expression level</b> <ul style="list-style-type: none"> <li>• Comparative transcriptomics</li> <li>• Disease biomarker</li> <li>• Homogenous systems</li> </ul>


Chaudhry et al. *scRNA-seq of the Cardiovascular System* (2019)

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

Thomas Conrad et al. (2020). *Maximizing transcription of nucleic acids with efficient T7 promoters*. *Communications Biology*, 3. 439.

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### Feature of LUTHOR

- High Quality Performance for **One Cell and Ultra-low input RNA**
- Direct Counting for **Gene Expression Quantification**
- One cell – One Comprehensive
- **THOR<sup>+</sup>** Amplification
- Direct RNA Amplification method


# Feature distribution of uniquely mapped reads for LUTHOR 3'mRNA-Seq. The majority of LUTHOR reads generated from 10 and 100pg UHRR, or 1 and 10 single HEK293T cells map to exonic sequences.

Method	Input	Detected Genes*
LUTHOR 3' mRNA-Seq	1 cell (HEK293T)	~11,000 – 15,000
SMART-Seq v4 3' DE	1 cell (K562)	~ 6,000 – 8,000 <sup>1</sup>
SMART-Seq Single Cell (WTS)	1 cell (GM22601)	~9,000 – 10,500 <sup>1</sup>
CEL-Seq 2 (C1 HT) (WTS)**	1 cell (mES)	~7,000 <sup>2</sup>
SMART-Seq 2 (WTS)**	1 cell (mES)	~8,000 – 9,000 <sup>2</sup>
SMART-Seq 3 (WTS)**	1 cell (HEK293FT)	~11,000 – 12,000 <sup>3</sup>

\*Obtained from manufacturers specifications [1] and from published benchmarking studies [2,3] at 1 M reads/sample, threshold of >1 Counts Per Million (CPM). Methods are shown for 3'-Seq and whole transcriptome sequencing (WTS). \*\* Academic protocol.

### "Single-Cell" RNA-Seq

Merging of Gene Expression Profiles of Hundreds of Cells

### LUTHOR RNA-Seq


One cell - One Comprehensive Expression Profile

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< Ebiogen's RNA-Seq Guide Map >

< 참고 문헌 >

1. 뇌종양에서의 단일 세포 RNA 시퀀싱을 이용한 전사체 연구동향. <http://www.ibric.org/myboard>. 2018, T37
2. Highly parallel genome-wide expression profiling of individual cells using nanoliter droplets. *Cell*\_2015, 161(5), 1202-1214
3. Droplet barcoding for single-cell transcriptomics applied to embryonic stem cells. *Cell*\_2015, 161(5), 1187-2101
4. Single-cell RNA sequencing 기술동향. <http://www.ibric.org/myboard>. 2018, T28
5. mRNA-Seq whole-transcriptome analysis of a single cell. *Nature Methods*\_2009,6,337-382
6. Maximizing transcription of nucleic acids with efficient T7 promoters. *Communications Biology*\_2020, 3, 439
7. TECH Note: Unprecedented sensitivity with the SMART-Seq Single Cell Kit. <https://www.takarabio.com>, 2020
8. Comparative analysis of Single-cell RNA Sequencing Methods. *Mol Cell*\_2017, 65, 631-643
9. Single-cell RNA counting at allele- and isoform-resolution using Smart-seq3. *Nat Biotechnol*\_2020, 38, 708-714
10. Single-Cell RNA Sequencing of the Cardiovascular System : New Looks for Old Diseases. *Front Cardiovasc Med*\_2019, 6, 173
11. Dissecting Cellular Heterogeneity Using Single-Cell RNA Sequencing. *Mol Cells*\_2019. 42(3), 189-199