
SINGLE CELL RNA-SEQ

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Sample requirement	Cell line, Primary cell, FACS sorted cell : 250,000 cells (Viability > 70%) Tissue ⁺ : Fresh Tissue 200 mg 이상 (in MACS buffer)
Library method	10X Genomics Next GEM-X 3' technology
Number of Target cell	3,000 ~ 20,000 cells
Turnaround Time	~ 4 weeks after Cell counting & QC
Data analysis	ExSCEA Report, UMAP, t-SNE, PCA, Clustering heatmap (Seurat/Loupe 기반 분석 제공 등)

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Sample requirement	Cell line, Primary cell, FACS sorted cell : 250,000 cells (Viability > 70%) Tissue ⁺ : Fresh Tissue 200 mg 이상 (in MACS buffer)
Library method	10X Genomics Next GEM-X 3' technology
Number of Target cell	5,000 cells
Turnaround Time	~ 4 weeks after Cell counting & QC
Data analysis	ExSCEA Report, UMAP, t-SNE, PCA, Clustering heatmap (Seurat/Loupe 기반 분석 제공 등)

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Sample requirement	Cell line, Primary cell, FACS sorted cell : 500,000 cells (Viability > 70%) Tissue ⁺ : Fresh Tissue 200 mg 이상 (in MACS buffer)
Library method	Illumina single cell 3' RNA-Seq preparation
Number of Target cell	10,000 cells / 20,000 cells
Turnaround Time	~ 4 weeks after Cell counting & QC
Data analysis	QC Report, Gene expression data, UMAP, t-SNE, Feature selection PCA, Clustering heatmap 등 (Dragen 기반 분석 제공)

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1. Hedlund E, Deng Q. Single-cell RNA sequencing: Technical advancements and biological applications. *Mol Aspects Med* 2018;**59**:36–46. <https://doi.org/10.1016/j.mam.2017.07.003>
2. Wang Z, Gerstein M, Snyder M. RNA-Seq: A revolutionary tool for transcriptomics. *Nat Rev Genet* 2009;**10**:57–63. <https://doi.org/10.1038/nrg2484>
3. Hwang B, Lee JH, Bang D. Single-cell RNA sequencing technologies and bioinformatics pipelines. *Exp Mol Med* 2018;**50**:1–14. <https://doi.org/10.1038/s12276-018-0071-8>
4. WOO, Hyunmin; EYUN, Seong-il. Applications and techniques of single-cell RNA sequencing across diverse species. *Briefings in bioinformatics*, 2025, 26.4: bbaf354. <https://doi.org/10.1093/bib/bbaf354>
5. NETLA, Vamsidhar Reddy, et al. A comparative analysis of single-cell transcriptomic technologies in plants and animals. *Current Plant Biology*, 2023, 35: 100289. <https://doi.org/10.1016/j.cpb.2023.100289>
6. DE SIMONE, Marco, et al. Comparative analysis of commercial single-cell RNA sequencing technologies. *BioRxiv*, 2024, 2024.06.18.599579. <https://doi.org/10.1101/2024.06.18.599579>
7. SWAMINATH, Sharmada; RUSSELL, Alistair B. The use of single-cell RNA-seq to study heterogeneity at varying levels of virus–host interactions. *PLoS Pathogens*, 2024, 20.1: e1011898. <https://doi.org/10.1371/journal.ppat.1011898>
8. CLARK, Iain C., et al. Microfluidics-free single-cell genomics with templated emulsification. *Nature Biotechnology*, 2023, 41.11: 1557–1566. <https://doi.org/10.1038/s41587-023-01685-z>
9. LE, Jiayuan, et al. Single-cell multi-omics in cancer immunotherapy: from tumor heterogeneity to personalized precision treatment. *Molecular Cancer*, 2025, 24.1: 221. <https://doi.org/10.1186/s12943-025-02426-3>
10. Ma L, et al. Tumor cell biodiversity drives microenvironmental reprogramming in liver cancer. *Cancer Cell*, 2019, 36:418–e430416. <https://doi.org/10.1016/j.ccell.2019.08.007>.
11. Liu Y, et al. Immune phenotypic linkage between colorectal cancer and liver metastasis. *Cancer Cell*, 2022, 40:424–e437425. <https://doi.org/10.1016/j.ccell.2022.02.013>.
12. Schumacher TN, Schreiber RD. Neoantigens in cancer immunotherapy. *Science*, 2015, 348:69–74. <https://doi.org/10.1126/science.aaa4971>.
13. Gubin MM, Artyomov MN, Mardis ER, Schreiber RD. Tumor neoantigens: Building a framework for personalized cancer immunotherapy. *J Clin Invest*, 2015, 125:3413–21. <https://doi.org/10.1172/jci80008>.
14. Koşaloğlu-Yalçın Z, et al. Predicting T cell recognition of MHC class I restricted neoepitopes. *Oncoimmunology*, 2018, 7:e1492508. <https://doi.org/10.1080/2162402x.2018.1492508>.
15. Cracking cancer: Using long-read RNA sequencing for cancer neoantigen discovery. 2023, <https://www.pacb.com/blog/cracking-cancer-using-long-read-rna-sequencing-for-cancer-neoantigen-discovery/>.
www.e-biogen.com

16. GIL DEL ALCAZAR, Carlos R., et al. Immune escape in breast cancer during in situ to invasive carcinoma transition. *Cancer discovery*, 2017, 7.10: 1098-1115. <https://doi.org/10.1158/2159-8290.CD-17-0222>
17. CUI, Xiaolu, et al. Single-cell and spatial transcriptomic analyses revealing tumor microenvironment remodeling after neoadjuvant chemoimmunotherapy in non-small cell lung cancer. *Molecular cancer*, 2025, 24.1: 111. <https://doi.org/10.1186/s12943-025-02287-w>
18. LI, Yuanfang, et al. CAF-macrophage crosstalk in tumour microenvironments governs the response to immune checkpoint blockade in gastric cancer peritoneal metastases. *Gut*, 2025, 74.3: 350-363. <https://doi.org/10.1136/gutjnl-2024-333617>
19. XU, Yanfen, et al. Multimodal single cell-resolved spatial proteomics reveal pancreatic tumor heterogeneity. *Nature Communications*, 2024, 15.1: 10100. <https://doi.org/10.1038/s41467-024-54438-0>
20. JI, Andrew L., et al. Multimodal analysis of composition and spatial architecture in human squamous cell carcinoma. *cell*, 2020, 182.2: 497-514. e22. <https://doi.org/10.1016/j.cell.2020.05.039>
21. ZHANG, Linan, et al. Single-cell analysis reveals the stromal dynamics and tumor-specific characteristics in the microenvironment of ovarian cancer. *Communications Biology*, 2024, 7.1: 20. <https://doi.org/10.1038/s42003-023-05733-x>
22. 이바이오젠 홈페이지, <https://www.e-biogen.com/service/single-cell-rna-seq>