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(A) Whole-Transcriptome RNA-Seq of FFPE Renal Cell Carcinoma Samples


그림 2. High concordance of total RNA-seq-based gene expression profiles between matched FF and FFPE samples [10]

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(A) Systematic benchmarking of imaging spatial transcriptomics platforms in FFPE tissues

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5.1 Ebiogen's Sample Submission Guide – FFPE Sample

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5.2 QuantSeq 3' mRNA-Seq

Sample requirement	>1 μ g total RNA (Min. >10ng total RNA)
Library method	QuantSeq 3' mRNA-Seq Library Prep Kit FWD
NGS run format	NextSeq, SE100
Data yield	>10M read/sample
Turnaround time	~2 weeks after RNA QC
Sample type	Tissue, Cell, Whole Blood(Paxgene), Plant, Fungi, FFPE , Sorted cell, etc.

5.3 Total RNA-Seq

Sample requirement	>2 μ g total RNA
Library method	RiboCop rRNA Depletion kit + CORALL RNA-seq Library Prep kit
NGS run format	NextSeq, PE150
Data yield	> 9 Gb/sample
Turnaround time	~3 weeks after RNA QC
Sample type	Tissue, Cell, Whole Blood(Paxgene), Plant, Fungi, FFPE , etc.

5.4 Visium HD Spatial Gene Expression

Species	Human, Mouse
Number of Gene	Human ~18,000 genes, Mouse ~19,000 genes
Sample Type	FFPE , Fresh Frozen, Fixed Frozen Tissue
Capture Size	6.5 mm x 6.5mm (1spot/slide)
Turnaround Time	~5 weeks after RNA QC

5.5 GeoMx Spatial Transcriptome + Protein Profiling

Service Name	GeoMx WTA* Service	GeoMx IPA Service	GeoMx DPA Service
Species	Human, Mouse	Human	Human
Number of Gene, protein	Human ~18,000 gene Mouse ~20,000 gene	~ 570 proteins	~ 1200 proteins
Data Yield	~100 reads/ μm^2	50 reads/ μm^2	50 reads/ μm^2
Sample Type	Tissue (FFPE , Tissue microarray, Fresh frozen tissue, Core needle biopsy etc.,)		
Capture Size	14.6 x 36.2 mm ²		
Turnaround time	~3 weeks after ROI selection		
Publication	NanoString GeoMx Publication		

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1. Conde et al., *The Utilization of Formalin Fixed-Paraffin-Embedded Specimens in High Throughput Genomic Studies*. (PLoS One. 2017)
2. Baker et al., *The impact of tissue handling on RNA integrity and transcriptome profiling*. (Nature Methods, 2010)
3. NanoString Technologies. Spatial biology sample types: FFPE, fresh frozen, and tissue microarray (TMA). Available from: <https://www.nanostring.com/>
4. Kokkat TJ, et al. *Archived Formalin-Fixed Paraffin-Embedded (FFPE) Blocks: A Valuable Resource for Clinical and Translational Research*. (PLoS ONE. 2013)
5. Siegel MB, et al. *RNA-seq from archival FFPE breast cancer samples: molecular pathway fidelity and novel discovery*. (BMC Medical Genomics. 2019)
6. Coudry et al., *Successful application of microarray and next-generation sequencing to formalin-fixed, paraffin-embedded tissue*. (Journal of Molecular Diagnostics, 2013)
7. Shabihkhani et al. *The procurement, storage, and quality assurance of frozen blood and tissue biospecimens in pathology, biorepository, and biobank settings*. (Clin Biochem. 2014)
8. Jacobsen et al., *Comparison of whole transcriptome sequencing of fresh, frozen, and formalin-fixed, paraffin-embedded cardiac tissue*. (Sci Rep. 2023)
9. Zhao et al., *Comparison of RNA-Seq by poly(A) capture, ribosomal RNA depletion, and DNA microarray for expression profiling* (BMC Genomics, 2014)
10. Li P, Conley A, Zhang H, Kim HL. *Whole-Transcriptome profiling of formalin-fixed, paraffin-embedded renal cell carcinoma by RNA-seq* (BMC Genomics, 2014)
11. 10x Genomics. *When spatial gene expression meets FFPE tissue blocks: A modern-day love story*. 10x Genomics Blog, May 5, 2021. Available at <https://www.10xgenomics.com/blog/when-spatial-gene-expression-meets-ffpe-tissue-blocks-a-modern-day-love-story>
12. Wang H, Huang R, Nelson J, et al. *Systematic benchmarking of imaging spatial transcriptomics platforms in FFPE tissues*. (Nature Communications, 2025)