Target 영역에 따른 DNA sequencing (WGS, WES, Targeted seq) 기술

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E-biogen's Service Information

WGS (Whole Genome Sequencing)

Sample requirement	>2ug gDNA
Library method	Illumina TruSeq Nano DNA Library kit
NGS run format	NovaSeq 6000, PE150
Turnaround time	~7 weeks after DNA QC
Data Analysis	SNP(Single Nucleotide Polymorphism), INDEL (Insertion & Deletion), C NV(Copy Number Variation), 개별 샘플에 대한 Structural Variation 분석 지원

WES (Whole Exome Sequencing)

Sample requirement	>2ug gDNA
Library method	Agilent Sure Select kit
NGS run format	NovaSeq 6000, PE100
Turnaround time	~6 weeks after DNA QC
Data Analysis	SNP(Single Nucleotide Polymorphism), INDEL (Insertion & Deletion), C NV(Copy Number Variation), 개별 샘플에 대한 Structural Variation 분석 지원

Targeted panel Sequencing

Sample requirement	>2ug gDNA
Library method	NEB Ultra DNA Library kit
NGS run format	NextSeq 550, PE150bp
Turnaround time	~6 weeks after DNA QC (kit 제작 기간 별도)
Data Analysis	ExPADA report, Peak Annotation, Motif Discovery, IGV, BedGraph 등

< 참고 문헌 >

1. Shendure, Jay, et al. "DNA sequencing at 40: past, present and future." *Nature* 550.7676 (2017): 345-353.
2. Richards, S., et al., *Standards and guidelines for the interpretation of sequence variants: a joint consensus recommendation of the American College of Medical Genetics and Genomics and the Association for Molecular Pathology. Genet Med, 2015. 17(5): p. 405-24*
3. Kosuri, Sriram, and George M. Church. "Large-scale de novo DNA synthesis: technologies and applications." *Nature methods* 11.5 (2014): 499-507.
4. NGS 기반+유전자검사(심화용), 식품의약품안전평가원, 2022
5. CELEMICS, End-to-end Customized Target Enrichment Panel, <https://www.celemics.com/products/customized-ngs-panel/>
6. Wang, Fang, et al. "Initial whole-genome sequencing and analysis of the host genetic contribution to COVID-19 severity and susceptibility." *Cell discovery* 6.1 (2020): 83.