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
그림.1 Body fluids for the identification of potential protein biomarkers for pancreatic cancer. [2]

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## Recent Trends in Biomarker Discovery Analysis using Body Fluids

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
그림2. Diagnostic opportunities of cell-free DNA and RNA within prenatal screening and disease diagnostics. [4]

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## Proteomics-Driven Body Fluid-Based Biomarkers

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
그림 3. Depicting the general method for proteomic analysis of cervical cancer derived from various body fluids. [1]

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Body Fluid Type	Methodology and Protocol	Cohort	Key Findings	Extra Comments (Merits/Demerits)	References
Plasma	2D-DIGE (silver staining); MS/MS (MALDI-TOF) to identify DEPs, and further validation by ELISA and statistical analysis by ANOVA.	Healthy: 40 individuals; CSCC and CIN patients: 80 individuals.	Cytokeratin 19 is upregulated in both the CIN 3 and CSCC IV conditions and tetranectin downregulated in CSCC.	Identification of DEPs along different stages of cervical cancer progression helps in understanding and prognosis of cancer.	[38]
Serum	Weak cation method, exchange chromatography fractionation in conjunction with MALDI-TOF spectrometry; liquid chromatography-electrospray ionization tandem mass spectrometry, and enzyme-linked immunosorbent assay (ELISA).	Healthy: 50 individuals; patients before surgery: 39; patients after surgery: 28.	The three peaks (m/z: 2435.63, 2575.3, and 2761.79 Da) may serve as predictive serum biomarkers for cervical cancer (CC).	Each patient group has obvious variation as the combined effect of age, stage, and tumor type reduces the power of marker detection.	[47]
Mucous	SELDI-TOF (surface-enhanced laser desorption and ionization-time of flight mass spectrometry).	Samples were collected from women attending urban hospital colposcopy clinics who were enrolled as a part of the study of cervical neoplasia. Samples were collected at the time of colposcopy by absorption into two Weck-Cel® sponges from 2-6 women matched for ages and races.	Annexin, tropomyosin, 14-3-3 sigma, calreticulin, and anterior gradient protein were identified.	The short sample size and inaccuracy of sample collecting techniques lowered the number of proteins discovered	[28]
Menstrual fluid	Genomic DNA was extracted from the menstrual blood collected on a napkin using a QIAamp DNA Mini Kit. Two rounds of PCR reaction using My11 and My09 primers for HPV detection. Fischer's exact test to examine the association between the distribution of genotypes or alleles for the TAP polymorphisms.	Control: 137 individuals; CIN3, CIN1, CIN2: 265 patients.	TAP1 1333V and TAP1 D637G were detected in the menstrual blood samples. The genotypes AA, AG, and GG were detected at each polymorphic site in the patients and the risk of developing high-grade cervical neoplasia was reduced for the AG and GG phenotypes as compared to the AA genotype. The risk of developing high-grade CIN was reduced in the patients that had a G allele than in those with an A allele.	The findings in the study have high specificity, sensitivity, and positive predictive value for the HPV virus and have received positive responses from over 5000 women.	[52]
Cervicovaginal fluid	Label-free quantification method based on LC-MS/MS method followed by ELISA. The PLS-DA model for further statistical analysis.	Development set—healthy: 10 individuals; LSIL: 10 individuals; HSIL: 10 individuals; cancer: 10 individuals Validation set—healthy: 14; LSIL: 8; HSIL: 6; cancer: 5.	ACTN4, VTN, ANXA1, ANXA2, CAP1, MUC5B and PKM2 from the 27 differentially expressed proteins have been indicated as promising biomarkers for cervical cancer.	The comparatively high number of samples gives better and more accurate results and reduces the chances of false biomarker discovery. The samples were also better classified into further four subgroups providing a comparison basis amongst the four groups.	[53]
Urine	Label-free quantification-UPLC-MS/MS analysis of pooled samples protein-protein interaction (STRING), pathway enrichment analysis, and molecular functions from KEGG and GO. Sensitivity as potential biomarkers tested by Western blotting and statistical analysis like logistic regression, ROC and AUC.	Healthy: 13 individuals; cervical cancer: 24 individuals.	Five Proteins with molecular weight >100 kDa were identified as potential biomarkers—LRG1, MMRN1 (upregulated), S100A1, CT344, SERPIN 33 (downregulated).	Rather than conventional gel-based MS analysis, non-gel based LPQ-MS analysis could aid in finding the low molecular weight potential biomarkers present in trace amounts in urine.	[55]

그림 4. Biomarkers identified from different body fluids of cervical cancer patients. [1]

## Cell free RNA Biomarkers Discovery from Body Fluids

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Molecule	Phenotype/test	Tissue(s)	Technology	References
miRNA	NAFLD and NAFLD-related disease	Liver biopsy, blood, rodent liver, serum, plasma	TaqMan microRNA array, NGS, RT-qPCR, biochemical markers, ultrasound	(3, 87-92)
	Steatosis	Serum, plasma	NMR, ultrasound, TaqMan microRNA array, NGS	(93)
	Obesity	Blood, MIR-34a KO mice	RT-qPCR, ELISA, biochemical markers	(94,95)
	T1D	Human and murine islets	Biochemical markers, RT-qPCR	(96)
	T2D	ZDF rats, blood, murine liver	RT-qPCR, NGS, LC-MS/MS	(97-99)
mRNA	Diabetic retinopathy	Serum	RT-qPCR	(100)
	Oral cancer	Saliva	Microarray, RT-qPCR	(101)
	Lung cancer	Serum and bronchial lavage	RT-qPCR	(102,103)
	Prostate cancer	Plasma, PBMC, serum	RT-qPCR	(31)
	Breast cancer	Serum	RT-qPCR, PCR-based fluorescence microsatellite assay	(102, 104,105)
	Renal cell carcinoma	Urine	RT-qPCR	(5)
	Bladder cancer	Urine	Microarray	(106)
	Rectal cancer	Plasma	RT-qPCR	(107)
	Fetal genetic abnormalities	Plasma, placental tissue	RT-qPCR, microarray	(108-111)
	Sex determination	Plasma	RT-qPCR	(112)
lncRNA	Pre-eclampsia	Plasma	RT-qPCR	(113)
	Diabetic retinopathy	Whole blood	RT-qPCR	(114,115)
	Diabetic consequences	Cell culture, murine model, blood	WB, ELISA, RT-qPCR, NGS, cloning, cell assays, microscopy	(116-118)
	Diabetic retinopathy	Cell culture, murine model, blood	ELISA, RT-qPCR, microarray, WB, cloning, microscopy	(119,120)

cfRNA, cell-free RNA; ELISA, enzyme-linked immunosorbent assay; KO, knock out; LC-MS/MS, liquid chromatography tandem mass spectrometry; NAFLD, non-alcoholic fatty liver disease; NGS, next-generation sequencing; NMR, nuclear magnetic resonance; PBMC, peripheral blood mononuclear cells; RT-qPCR, real-time quantitative polymerase chain reaction; T1D, type 1 diabetes; T2D, type 2 diabetes, WB, Western blot; ZDF rat, Zucker diabetic fatty rat.

그림5. Examples of studies on circulating cell-free RNA as potential diagnostic tool. [4]

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### E-biogen's Body Fluids Analysis Service

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Service	Antibody array		*LC-MS Proteomics , PTM	**Small RNA-Seq, Exosomal small RNA-Seq	**Arraystar microarray
	RayBiotech	Full moon BioSystems			
Application	Protein profiling	Protein, phosphorylation profiling	Protein , PTM (glycosylation, phosphorylation) profiling	miRNA profiling	LncRNA, circular , piRNA, Super Enhancer LncRNA profiling
Sample Type	Serum/plasma, Exosome, other body fluids	Serum/plasma	Serum/plasma, Exosome CSF, urine, saliva, other body fluids	Whole blood, Serum/plasma , Exosome, Urine	Whole blood, Serum/plasma, Exosome
Experiment method	High-throughput ELISA Antibody based Array		Agilent 1260 & 1290 System Thermo Q-Exactive MS	NextSeq 500/550, SE 75bp	RNA microarray designed specific probes on chip
Data type	Semi-quantitative/ quantitative	Semi-quantitative	Semi-quantitative	Semi-quantitative	Semi-quantitative
Turnaround time	~ 4 weeks after Protein QC		~ 6 weeks after Protein QC	~ 4 weeks after RNA QC	~ 6 weeks after RNA QC
Link	<a href="#">Arrays   RayBiotech</a>	<a href="#">Arrays - Full Moon BioSystems</a>	<a href="#">ebiogen Inc. (e-biogen.com)</a>	<a href="#">ebiogen Inc. (e-biogen.com)</a>	<a href="#">Featured ncRNA Microarrays - piRNA Microarray - Super-enhancer LncRNA Array</a>

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< 참고문헌 >

1. Insights on Proteomics-Driven Body Fluid-Based Biomarkers of Cervical Cancer, *Proteomes*, 2022, 10, 13.
2. Proteomic biomarkers in body fluids associated with pancreatic cancer, *Oncotarget*, 2018, Vol. 9, (No. 23), pp: 16573-16587
3. High-throughput miRNA sequencing and identification of biomarkers for forensically relevant biological fluids, *Electrophoresis*, 2016, Volume 37, Issue 21, Pages 2780-2788
4. Cell-free DNA and RNA—measurement and applications in clinical diagnostics with focus on metabolic disorders, *Physiol Genomics* 53: 33–46, 2021.
5. Human body fluid proteome analysis, *Proteomics*. Author manuscript; available in PMC 2010 July 26. NIH-PA
6. Review of Liquid Chromatography-Mass Spectrometry-Based Proteomic Analyses of Body Fluids to Diagnose Infectious Diseases, *Int. J. Mol. Sci.* 2022, 23, 2187.
7. Body fluid biomarkers in multiple sclerosis, *Lancet Neurol* 2014; 13: 113–26
8. Current challenges and best practices for cell-free long RNA biomarker discovery, Cabús et al. *Biomarker Research* (2022) 10:62
9. Exosome-derived noncoding RNAs in gastric cancer: functions and clinical applications, Tang et al. *Mol Cancer* (2021) 20:99
10. Extracellular miRNAs: From Biomarkers to Mediators of Physiology and Disease, *Cell Metab.* 2019 October 01; 30(4): 656–673
11. MicroRNAs in body fluids—the mix of hormones and biomarkers, *Nat Rev Clin Oncol.* ; 8(8): 467–477.
12. Body Fluid Collection And Diagnostics Market Size, Share & Trends Analysis Report By Sample Type (Blood, Saliva, Urine, Cerebrospinal Fluid), By Products, By Technology, By Application, By Region, And Segment Forecasts, 2023 – 2030 Report ID: GVR-4-68040-019-8
13. HBFP: a new repository for human body fluid proteome, *Database*, Volume 2021, 2021, baab065