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
그림 4. The tissue morphology was delineated by the immunofluorescence detection of PanCK (epithelial cytokeratin, green), CD45+ (immune cells, magenta) proteins and SYTO13 (nuclei, blue) (B).

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
그림 5. Regions of interest (ROIs) were selected for molecular profiling analysis from the slides of a pedunculated tubular adenoma with low-grade dysplasia as a control (A) and pedunculated T1 colorectal cancers (B).

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
그림 6. A) GeoMx images with fluorescent morphology markers for CD45 (magenta), PanCK (green) and nuclei (blue). Samples were segmented on the GeoMx DSP for PanCK positive (green) and PanCK negative (red) areas of interest (AOI).

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	GeoMx Spatial WTA Service	GeoMx Spatial CTA Service	GeoMx Spatial IPA Service
Species	Human, Mouse	Human	Human
Target	Human 18,000개, Mouse 23,000개	Human 1,800개	Human 580개
Data Yield	12AOI 기준, 100reads/um ²	12AOI 기준, 30reads/um ²	12AOI 기준, 200reads/um ²
Sample type	Tissue (FFPE, Tissue Microarray, Fresh frozen, Core needle biopsy)		
Sample size	14.1 X 35.1mm		
Turnaround time	3 weeks after ROI selection		

*WTA (Whole Transcriptome Atlas) / *CTA (Cancer Transcriptome Atlas)/ *IPA(Immuno-Oncology Proteome Atlas)

<참고문헌>

- Rosenbloom, Alyssa, et al. "120 Human immuno-oncology proteome atlas enables a holistic proteomic approach for spatial biology discovery in FFPE tissues." (2023).
- Church, Sarah E., et al. "Translational and clinical applications of the GeoMx digital spatial profiling platform." *Diagnostic Molecular Pathology*. Academic Press, 2024. 767-783.
- Janesick, Amanda, et al. "High resolution mapping of the tumor microenvironment using integrated single-cell, spatial and in situ analysis." *Nature Communications* 14.1 (2023): 8353.
- <https://nanostring.com/products/geomx-digital-spatial-profiler/geomx-protein-assays/io-proteome-atlas/>
- Shamis, Suad Ahmad Khalifa, et al. "Spatial transcriptomic analysis of tumour with high and low CAIX expression in TNBC tissue samples using GeoMx™ RNA assay." *Histology and Histopathology: Cellular and Molecular Biology* (2023).
- Frei, Nicola F., et al. "Feasibility Study Utilizing NanoString's Digital Spatial Profiling (DSP) Technology for Characterizing the Immune Microenvironment in Barrett's Esophagus Formalin-Fixed Paraffin-Embedded Tissues." *Cancers* 15.24 (2023): 5895.
- <https://nanostring.com/products/geomx-digital-spatial-profiler/geomx-dsp-overview/>
- https://www.e-biogen.com/board_tech.php?act=view&aid=1333&page=1
- https://www.e-biogen.com/board_tech.php?act=view&aid=888&page=1
- Kwon, Mi Jung, et al. "Potential Molecular Markers Related to Lymph Node Metastasis and Stalk Resection Margins in Pedunculated T1 Colorectal Cancers Using Digital Spatial Profiling: A Pilot Study with a Small Case Series." *International Journal of Molecular Sciences* 25.2 (2024): 1103.
- Chumsri, Saranya, et al. "207 The spatially-resolved proteomic atlas of the tumor immune architecture of HER2-positive breast cancer and response to neoadjuvant anti-HER2 therapy-based chemotherapy." (2023).
- Tanjak, Pariyada, et al. "The KRAS-mutant consensus molecular subtype 3 reveals an immunosuppressive tumor microenvironment in colorectal cancer." *Cancers* 15.4 (2023): 1098.