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EM-Seq (Enzymatic Methyl-Sequencing)

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(A) When nucleosomes are spaced closely together, transcription factors cannot bind and gene expression is turned off.

(B) When nucleosomes are spaced far apart, transcription factors can bind, allowing gene expression to occur.[3]

ChIP-Seq (Chromatin Immunoprecipitation Sequencing)

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ATAC-Seq (Assay for Transposase-Accessible Chromatin using Sequencing)

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Integrative Analysis of ATAC-Seq and RNA-Seq Data

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
그림 8. Differential Expression and Functional Enrichment of Candidate DEGs. [18]

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Integrative Analysis of ChIP-Seq and RNA-Seq Data

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
그림 9. PML/RAR α exerts both repressive and activating functions through direct binding. [19]

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Ebiogen's Epigenome and Transcriptome Sequencing Service

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Enzymatic Methyl-Seq(EM-Seq) Service

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Sample requirement	>2ug gDNA
Sample Type	gDNA
Library method	NEBNext® Enzymatic Methyl-seq Kit
NGS run format	Illumina NovaSeq6000, PE150bp
Data yield	~120G/sample
Turnaround time	~6 weeks after DNA QC

Chromatin Profiling (ChIP-Seq) Service

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Sample requirement	>20ul IP-DNA
Sample Type	IP-DNA
Library method	NEB Ultra DNA Library kit
NGS run format	Illumina NovaSeq6000, PE100bp
Data yield	~4G/sample
Turnaround time	~3 weeks after Library QC

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Sample requirement	> ~1x10 ⁶ cells in 1ml media (Cell viability >80%)
Sample Type	Cell
Library method	ACTIVE MOTIF® ATAC-Seq Kit
NGS run format	Illumina NovaSeq6000, PE100bp
Data yield	~6G/ sample
Turnaround time	~5 weeks after Cell QC

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	Total RNA-Seq	mRNA-Seq
Sample requirement	>2ug total RNA	>2ug total RNA
Library method	RiboCop rRNA Depletion kit	Poly(A) RNA Selection kit

	+ CORALL RNA-seq Library Prep kit	+ CORALL RNA-seq Library Prep kit
NGS run format	HiSeq or NovaSeq, PE100	
Data yield	>6Gb/sample	>4Gb/sample
Turnaround time	~4 weeks after RNA QC	
Sample type	Tissue, Cell, Whole Blood(Paxgene), Plant, Fungi, Bacteria, etc.	

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Sample requirement	>10ng total RNA
Sample Type	Tissue, Cell, Whole Blood(Paxgene), Plant, Fungi, FFPE, Sorted cell, etc.
Library method	QuantSeq 3' mRNA-Seq Library Prep Kit FWD
NGS run format	NextSeq 500/550, SE 75bp
Data yield	>10M read/sample
Turnaround time	~3 weeks after RNA QC

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ExPADA(Excel based Peak Annotation and Differential enrichment Analysis)

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< 참고 문헌 >

1. Illumina. (2022). *Multimomics eBook*
2. Wikipedia. (2023). *Regulation of gene expression*.
3. RWU Pressbooks. (2023). *Chapter 17. Regulation of Gene Expression – Introduction to Molecular and Cell Biology*.
4. Alokail, M.S., & Alenad, A.M. (2015). *DNA Methylation*. In M. Gunduz (Ed.), *A Concise Review of Molecular Pathology of Breast Cancer*. IntechOpen.
5. Cui, G., Dong, Q., Gai, K., & Qi, S. (2022). *Chromatin Dynamics: Chromatin Remodeler, Epigenetic Modification and Diseases*. In T. Huang (Ed.), *Epigenetics - Regulation and New Perspectives*. IntechOpen
6. International Human Genome Sequencing Consortium. (2001). *Initial sequencing and analysis of the human genome*. *Nature*, 409, 860–921.
7. Schübeler D. 2015. *Function and information content of DNA methylation*. *Nature* 517: 321–326. 10.1038/nature14192
8. Smith ZD, Meissner A. (2013). *DNA methylation: roles in mammalian development*. *Nat Rev Genet* 14: 204–220.
9. Vaisvila R, Ponnaluri VKC, Sun Z, et al. (2021). *Enzymatic methyl sequencing detects DNA methylation at single-base resolution from picograms of DNA*. *Genome Res* 31: 1280–1289.
10. Skvortsova, K., Zotenko, E., & Luu, P.-L. (2017). *Comprehensive evaluation of genome-wide 5-hydroxymethylcytosine profiling approaches in human DNA*. *Epigenetics & Chromatin*, 10, Article 16.
11. New England Biolabs. (n.d.). *NEBNext Enzymatic Methyl-seq (EM-seq) Technical Note*
12. Nakato, R., & Shirahige, K. (2017). *Recent advances in ChIP-seq analysis: from quality management to whole-genome annotation*. *Briefings in Bioinformatics*, 18(2), 279-290.
13. Sundaram, A. Y. M., Hughes, T., Biondi, S., et al. (2016). *A comparative study of ChIP-seq sequencing library preparation methods*. *BMC Genomics*, 17, Article 816.
14. Cell Signaling Technology. [ChIP Protocol Steps]. Retrieved from [<https://www.cellsignal.com/applications/chip-and-chip-seq/regulation-expression-in-cell-and-tissue>]
15. 이바이오젠. (2022). [기술노트] *Epigenetics (후성유전학) 개요 및 기술*.
16. Itai, Y., Rappoport, N., & Shamir, R. (2023). *Integration of gene expression and DNA methylation data across different experiments*. *Nucleic Acids Research*, 51(15), 7762–7776.
17. Hermani, S., Yarmolinsky, J., & Millstein, J. (2018). *Making sense of the epigenome using data integration approaches*. *Frontiers in Genetics*.
18. Chen, N., Liu, Q., Ma, X., et al. (2021). *Integrative Analysis of ATAC-Seq and Whole-Transcriptome Sequencing in the Trilog of Gastric Carcinogenesis: From Normal Mucosa to Precancerous Lesions and to Gastric Cancer*.
19. Tan, Y., Wang, X., & Song, H. (2021). *A PML/RARa direct target atlas redefines transcriptional deregulation in acute promyelocytic leukemia*. *Blood*, 137(11), 1503-1516.