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Shotgun Metagenome Seq.

Sample Requirement	>5µg gDNA, >50ng/µl gDNA, >4 DIN
Library Method	NEB Ultra DNA Library kit
NGS run format	NovaSeq, PE150bp
Data yield	~10Gb/sample
Turnaround time	~ 4 weeks after Library QC
Sample type	gDNA (prep 서비스도 제공)

Metagenome (16S V3-V4 Amplicon Seq.)

Sample Requirement	>0.1 ng/µl gDNA
Library Method	16S rRNA Amplicon Library(16S/ITS) for Illumina
NGS run format	MiSeq, PE300bp
Data yield	~50,000 reads/sample
Turnaround time	~ 4 weeks after Library QC
Sample type	gDNA (prep 서비스도 제공)

Metagenome (16S rRNA Full length Seq.)

Sample Requirement	>0.1 ng/μl gDNA
Library Method	SMRTbell prep kit for PacBio
NGS run format	PacBio Revio / Sequel I & Sequel II / RS II
Data yield	~50,000 reads/sample
Turnaround time	~ 4 weeks after Library QC
Sample type	gDNA (prep 서비스도 제공)

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Untargeted Metabolomics

Sample Type	Cell pellet, Tissue, Serum/Plasma, Bacteria, Fecal, Plant, etc.
Analysis Instrument	LC: Agilent 1290 MS: Thermo Q-Exactive MS (High Resolution MS)
Software	Compound Discoverer, MetaboAnalyst
Turnaround time	~ 7 weeks after sample arrival

< 참고 문헌 >

1. <https://en.wikipedia.org/wiki/Metagenomics>
2. Thomas T, Gilbert J, Meyer F. *Metagenomics - a guide from sampling to data analysis*. *Microb Inform Exp*. 2012 Feb 9;2(1):3. doi: 10.1186/2042-5783-2-3. PMID: 22587947; PMCID: PMC3351745.
3. Joynson, R., Pritchard, L., Osemwekha, E., & Ferry, N. (2017). *Metagenomic Analysis of the Gut Microbiome of the Common Black Slug Arion ater in Search of Novel Lignocellulose Degrading Enzymes*. *Frontiers in Microbiology*, 8. doi:10.3389/fmicb.2017.02181
4. The Integrative HMP (iHMP) Research Network Consortium. *The Integrative Human Microbiome Project*. *Nature* 569, 641–648 (2019). <https://doi.org/10.1038/s41586-019-1238-8>
5. Ranjan, R., Rani, A., Metwally, A., McGee, H. S., & Perkins, D. L. (2016). *Analysis of the microbiome: Advantages of whole genome shotgun versus 16S amplicon sequencing*. *Biochemical and biophysical research communications*, 469(4), 967-977.
6. https://www.e-biogen.com/support/tech-info?tpf=board/view&board_code=2&code=2757
7. Sharpton, T. J. (2014). *An introduction to the analysis of shotgun metagenomic data*. *Frontiers in plant science*, 5, 209.
8. Quince, C., Walker, A. W., Simpson, J. T., Loman, N. J., & Segata, N. (2017). *Shotgun metagenomics, from sampling to analysis*. *Nature biotechnology*, 35(9), 833-844.
9. https://astrobiomike.github.io/misc/amplicon_and_metagen
10. <https://www.e-biogen.com/service>
11. Bolger, A.M., Lohse, M. & Usadel, B. *Trimmomatic: a flexible trimmer for Illumina sequence data*. *Bioinformatics* 30, 2114–2120 (2014).
12. Strous, M., Kraft, B., Bisdorf, R., & Tegetmeyer, H. E. (2012). *The binning of metagenomic contigs for microbial physiology of mixed cultures*. *Frontiers in microbiology*, 3, 410.
13. Abubucker, S. et al. *Metabolic reconstruction for metagenomic data and its application to the human microbiome*. *PLoS Comput. Biol.* 8, e1002358 (2012).
14. Aguiar-Pulido, V., Huang, W., Suarez-Ulloa, V., Cickovski, T., Mathee, K., & Narasimhan, G. (2016). *Metagenomics, metatranscriptomics, and metabolomics approaches for microbiome analysis: supplementary issue: bioinformatics methods and applications for big metagenomics data*. *Evolutionary bioinformatics*, 12, EBO-S36436.
15. Lai, W. T., Deng, W. F., Xu, S. X., Zhao, J., Xu, D., Liu, Y. H., ... & Rong, H. (2021). *Shotgun metagenomics reveals both taxonomic and tryptophan pathway differences of gut microbiota in major depressive disorder patients*. *Psychological medicine*, 51(1), 90-101.
16. Enebe, M. C., & Babalola, O. O. (2020). *Effects of inorganic and organic treatments on the microbial community of maize rhizosphere by a shotgun metagenomics approach*. *Annals of Microbiology*, 70, 1-10.