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### Targeted vs Untargeted Metabolomics

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표 1. Targeted Metabolomics Untargeted Metabolomics 의 장단점 [6]

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
그림 3. Findability, Accessibility, Interoperability, Reusability, FAIR[11]

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
그림 5 Plasma metabolites dysregulations in patients with metastatic pNETs vs patients with non-metastatic pNETs[12]


그림 6 Differential plasma metabolites in pNET patients with and without MEN1 mutation[12]

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
그림 7. EPIC (European Prospective Investigation into Cancer and Nutrition) Consortium Network Overview[3]

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
그림 11. Notable populationbased studies analyzed by biocrates technology for metabolic phenotyping[3]

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## Metabolite coverage

MxP® Quant 500 kit		
Alkaloids (1)	Hormones (4)	Hexosylceramides (19)
Amine oxides (1)	Indoles and derivatives (4)	Dihexosylceramides (9)
Amino acids (20)	Nucleobases and related (2)	Trihexosylceramides (6)
Amino acid related (30)	Vitamins and cofactors (1)	Cholesteryl esters (22)
Bile acids (14)	Acylcarnitines (40)	Diglycerides (44)
Biogenic amines (9)	Phosphatidylcholines (76)	Triglycerides (242)
Carbohydrates and related (1)	Lysophosphatidylcholines (14)	
Carboxylic acids (7)	Sphingomyelins (15)	
Cresols (1)	Ceramides (28)	
Fatty acids (12)	Dihydroceramides (8)	

<b>MxP® Quant 500 XL kit</b>		
Alkaloids (1)	Vitamins and cofactors (1)	Sphinganine and sphingosine (8)
Amine oxides (1)	Phosphatidic acids (41)	Sphinganine and sphingosine phosphates (8)
Amino acids (20)	Lysophosphatidic acids (8)	Sphingomyelins (15)
Amino acid related (30)	Phosphatidylcholines (76)	Ceramides (29)
Bile acids (14)	Lysophosphatidylcholines (14)	Dihydroceramides (8)
Biogenic amines (9)	Phosphatidylethanolamines (95)	Hexosylceramides (19)
Carbohydrates and related (1)	Lysophosphatidylethanolamines (43)	Dihexosylceramides (9)
Carboxylic acids (7)	Phosphatidylglycerols (64)	Trihexosylceramides (6)
Cresols (1)	Lysophosphatidylglycerols (10)	Acylcarnitines (40)
Fatty acids Free/non-covalently bound (12)	Phosphatidylinositols (53)	Cholesteryl esters (22)
Hormones and related (4)	Lysophosphatidylinositols (16)	Monoglycerides (12)
Indoles and derivatives (4)	Phosphatidylserines (18)	Diglycerides (44)
Nucleobases and related (2)	Lysophosphatidylserines (12)	Triglycerides (242)

<b>AbsoluteIDQ p180 kit</b>		
Amino acids (21)	Acylcarnitines (40)	Sphingolipids (15)
Biogenic amines (21)	Lysophosphatidylcholines (14)	
Hexose (1)	Phosphatidylcholines (76)	

<b>AbsoluteIDQ Bile Acids kit</b>	
10 primary bile acids	produced by the liver (CA, GCA, TCA, CDCA, GCDCA, TCDCA, $\alpha$ -MCA, $\beta$ -MCA, $\omega$ -MCA, TMCA( $\alpha$ + $\beta$ ))
10 secondary bile acids	produced by gut bacteria from primary bile acids (DCA, GDCA, TDCA, LCA, GLCA, TLCA, UDCA, GUDCA, TUDCA, HDCA)

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<b>Sample Type</b>	20µL of Cell lysate, Serum/Plasma, Media, Homogenate (tissue, fecal, etc.)
<b>Kit Types</b>	MxP® Quant 500 kit MxP® Quant 500 XL kit AbsoluteIDQ p180 kit AbsoluteIDQ Bile Acids kit
<b>Analysis Instrument</b>	Thermo Orbitrap Eclipse Tribrid
<b>Software</b>	WebIDQ
<b>Turnaround time</b>	~ 7 weeks after Sample Arrival

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<b>Sample Type</b>	Cell pellet, Tissue, Serum/Plasma, Exosome pellet, Bacteria pellet, Plant, Fecal, etc.
<b>Analysis Instrument (HPLC &amp; UPLC)</b>	Agilent 1290 LC Thermo Vanquish LC
<b>Analysis Instrument (Mass Spectrometry)</b>	Thermo Q-Exactive MS (High Resolution MS) Thermo TSQ Altis (MRM)
<b>Software</b>	Compound Discoverer, MetaboAnalyst
<b>Turnaround time</b>	~ 7 weeks after Sample Arrival

< 참고 문헌 >

1. Qiu S, Cai Y, Yao H, Lin C, Xie Y, et al. Small molecule metabolites: discovery of biomarkers and therapeutic targets. Signal Transduct Target Ther. 2023.
2. Metabolon. Your Guide to Metabolomics Chapter 1—Metabolomics, Metabolites, and the Metabolome.
3. Biocrates life sciences ag. Whitepaper. FAIR compliant metabolomics profiling of population-based studies.
4. biocrates life sciences ag. Why metabolomics? <https://biocrates.com/why-metabolomics/>
5. Ebiogen. Technical Note, Jul. 2023. 다양한 분야에서 응용가능한 LC-MS/MS Metabolomics 분석 기술
6. Metabolon. Blog. Targeted vs Untargeted Metabolomics
7. Thomas S, Targeted vs. Untargeted Metabolomics. Ometa labs.
8. Chen, Y., Wang, B., Zhao, Y., et al. Metabolomic machine learning predictor for diagnosis and prognosis of gastric cancer. nature communications. 2024.
9. López-Hernández, Y., Monárrez-Espino, J., Herrera-van Oostdam, A., et al. Targeted metabolomics identifies high performing diagnostic and prognostic biomarkers for COVID-19. scientific reports. 2021.
10. Saigusa D, Hishinuma E, Matsukawa N, et al. Comparison of Kit-Based Metabolomics with Other Methodologies in a Large Cohort, towards Establishing Reference Values. Metabolites. 2021.
11. National Library of Medicine. Definition of FAIR Data. <https://www.nlm.nih.gov/oet/ed/cde/tutorial/02-200.html>
12. Jannin A, Dabo-Niang S, Do Cao C, et al. Identification of metabolite biomarkers for pancreatic neuroendocrine tumours using a metabolomic approach. European Journal of Endocrinology. 2025
13. Langsdorf, M., Ömer, G., Dearth, S., et al. Standardized quantitative metabolomics using the biocrates MxP® Quant 500 XL kit across mass spectrometer platforms. Biocrates. Application Note v1-2023.