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
그림 1. Single-cell RNA-sequencing, single-nuclei RNA-sequencing and spatial transcriptomics studies reveal cellular and molecular heterogeneity in human neurological disorders [2]

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
그림 2. Experimental pipeline of a droplet-based single-cell RNA sequencing in a nutshell. [2]

	Single cell RNA-Seq	Single nucleus RNA-Seq
Pros	• Captures cytoplasmic and nuclear transcripts • Enhanced capture of immune cells for disease studies • Large number of datasets have been generated	• Can be used on archived samples (fixed or frozen) • Improved capture of hard to dissociate or fragile cells
Cons	• Tissue needs to be fresh • Dissociation induces stress response which can impact data interpretation • Biases toward cells that survive dissociation process	• Depletion of immune cells • Capture of mostly unspliced pre-mRNA

표 1. Comparing information from scRNAseq vs snRNAseq [3]

Single-cell technology in neuroscience

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
그림 3. High-dimensional single-cell RNA sequencing data can be visualized by using dimensionality reduction algorithms to reveal cell clusters [2]

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
Method	Considerations	Publication/Protocol
Density gradient	Requires large sample input, as significant cell loss is expected	• Nguyen, H. X., Beck, K. D., & Anderson, A. J. (2011). • Nikodemova, M., & Watters, J. J. (2012). • Del-Aguila et al (2019).
Myelin removal kit	Requires large sample input	• Isolation of Nuclei for Single Cell RNA Sequencing & Tissues for Single Cell RNA Sequencing • Miltenyi Myelin Removal Beads
FACS	Useful when trying to enrich for specific cell types and cleanup at the same time Requires large sample input	• Pan, J., & Wan, J. (2020)

표 2. Myelin Removal method


Spatial transcriptomics technology

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a Major events in evolution of prequel techniques

b Major prequel WM ISH atlases

c Major events in evolution of current-era techniques

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
그림 8. A side-by-side comparison of Visium v2 data (left) and Visium HD data (right) in FFPE human colorectal cancer, demonstrating the enhanced discovery power of whole transcriptome spatial gene expression at single cell-scale resolution.[8]

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Spatial transcriptomics technology in neuroscience

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Service Info.



Sample requirement+	-Fresh Tissue 200 mg++ - Cell line : $\sim 1 \times 10^6$ cells in 14ml media(Cell viability >90%) - Primary cell or FACS sorted cell : 250,000 cells in 500 ul media(Viability >70%)
Library method	10X Genomics Next GEM technology
NGS run format	Nova-Seq 6000, PE100bp
Data yield+	$\sim 60,000$ reads(30,000 Paired Reads) per cell
Number of Target cell+	$\sim 5,000$ cells
Turnaround time	~ 4 weeks after Cell counting & QC
Sample type	Cell, Nuclei, Tissue

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


Service Info.

Service Name	GeoMx Spatial WTA* Service	GeoMx Spatial CTA* Service	GeoMx Spatial IPA* Service
Species	Human, Mouse	Human	Human
Number of Gene	Human ~18,000ea / Mouse ~22,000ea	~1,800ea	~570ea
Data yield	~100 reads/ μm^2	~30 reads/ μm^2	~200 reads/ μm^2
Sample type	Tissue (FFPE, Tissue microarray, Fresh frozen tissue, Core needle biopsy etc.,)		
Sample size	14.6 x 36.2 mm		
Turnaround time	~3 weeks after ROI selection		
Publication	NanoString GeoMx Publication		

*WTA (Whole Transcriptome Atlas) / *CTA (Cancer Transcriptome Atlas) / *IPA (Immuno-Oncology Proteome Atlas)

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<참고 문헌>

1. Jung, N. & Kim, T.-K. Spatial transcriptomics in neuroscience. *Exp. Mol. Med.* 55, 2105–2115 (2023).
2. Piwecka, M., Rajewsky, N. & Rybak-Wolf, A. Single-cell and spatial transcriptomics: deciphering brain complexity in health and disease. *Nat. Rev. Neurol.* 19, 346–362 (2023).
3. [Comparative Analysis of Single-cell and Single-nucleus RNA Sequencing](#)
4. Jill M Williamson & David A Lyons Myelin Dynamics Throughout Life: An Ever-Changing Landscape? *Front Cell Neurosci.* 2018 Nov 19:12:424.
5. Moses L, Pachter L. Museum of spatial transcriptomics. *Nat Methods.* 2022;19(5):534-546. doi:10.1038/s41592-022-01409-2
6. <https://nanosttring.com/products/geomx-digital-spatial-profiler/geomx-dsp-overview/>
7. [Visium HD: Spatially Resolved, Single Cell Scale Gene Expression Mapping of FFPE Tissue](#)
8. [Your introduction to Visium HD: Spatial biology in high definition](#), 10X genomics Blog